

CR500D/CR1000D

Digital Wheelchair Scale



Featuring built-in handles, wheelstop, and wheels, the Roll-A-Weigh CR500D/CR1000D can be moved wherever you need it. Whether it's multiple weighing locations at a single facility or a stationary location, the Roll-A-Weigh gives you portability and ease of use and takes six "C" batteries or optional AC adapter.

Capacity:**CR500D**

500 lb x 0.2 lb/225 kg x 0.1 kg

CR1000D

1,000 lb x 0.5 lb/450 kg x 0.2 kg

Overall Dimensions:

32.5" W x 42.5" D x 3.5" H/83 cm W x 108 cm D x 9 cm H

Platform Size:

30" W x 32.5" D/76 cm W x 83 cm D



The perfect scale for weighing patients in manual or power wheelchairs, dialysis patients, or walk-ons. The lightweight Roll-A-Weigh features built-in wheels for easy portability.

Detecto's Roll-A-Weigh Wheelchair Scale features a weight display which can be mounted on the scale, or remotely on a wall or desk for convenient viewing.