



5079-325
Harvey™ DLX Double Head Stethoscope



5079-321
Harvey™ DLX Triple Head Stethoscope

Harvey™ DLX Stethoscopes

Developed for precision and acoustic clarity.

- > Optimized for listening to high-frequency heart sounds including murmurs, clicks, and ejection sounds
- > Bell ideally sized to bridge intercostal space between ribs for excellent clarity of low- and medium-frequency sounds
- > Triple-head models also have corrugated diaphragm for general overview, or for detecting low-frequency sounds, rumbles, and gallops
- > Brass chestpiece warms easily and provides superior acoustic response
- > Adjustable binaurals allow you to custom fit your stethoscope for optimal performance and comfort
- > Latex-free

	5079-321	Harvey™ DLX, Triple Head, Black
	5079-322	Harvey™ DLX, Triple Head, Burgundy
	5079-323	Harvey™ DLX, Triple Head, Navy
	5079-324	Harvey™ DLX, Triple Head, Forest Green
	5079-325	Harvey™ DLX, Double Head, Black
	5079-326	Harvey™ DLX, Double Head, Burgundy
	5079-327	Harvey™ DLX, Double Head, Navy
	5079-328	Harvey™ DLX, Double Head, Forest Green



1-866-588-3888 www.cwimedical.com