

ID field for electronic labeling of thermometer



Battery life indicator

User-selectable icon and/or words for identification of thermometer mode

Interchangeable probe with push-button probe cover ejection

Removable probe well for easy cleaning

Convenient storage for 25 probe covers

C/F button for toggling between Fahrenheit and Celsius

Large, easy-to-read, backlit LCD

Mode button for selecting oral, pediatric axillary, or adult axillary

Button for timing pulse and respiration

Ergonomic design with rubberized grips

Prominent recall button for display of last temperature taken

Storage space for 50 additional probe covers in the wall mount (Wall mount standard with the 692 model)

Interchangeable rectal probe and well reduce risk of cross-contamination



## SureTemp® Plus 692 specifications

Calibration Accuracy ..... ±0.2 °F (±0.1 °C) (Monitor Mode)

Temperature Predict Time:

Oral ..... Approx. 4-6 secs

Adult Axillary (18 years and older) ..... Approx. 15 secs

Pediatric Axillary (17 years and younger) ..... Approx. 10 secs

Rectal ..... Approx. 10 secs

Unit Dimensions ..... 8.46" x 3.18" x 2.43"  
(215 mm x 81 mm x 62 mm)

Weight ..... 12.6 ounces (357 grams)

Ambient Temperature Range ..... 50.0 °F (10 °C) - 104 °F (40.0 °C)

Patient Temperature Range ..... 80.0 °F (26.7 °C) - 110 °F (43.3 °C)

Power Source ..... (3) 1.5 Vdc AA batteries

Battery Life ..... Approx. 5000 temperatures

Warranty ..... 3-year limited instrument,  
1-year limited probe, and  
90-day limited probe well



**CWI Medical**

Medical Supplies & Health Care Products

1-866-588-3888  
[www.cwimedical.com](http://www.cwimedical.com)